

CBS-300 Electrolyte Analyzer (K^+ Na^+ Cl^-)

WWW.INSIGNIABIMEDICAL.COM



Products Characteristics

- Advanced Multi-Biosensor Technology(MBT), Free of Maintenance.
- Close Reagent Pack, Easy to Change
- Durable Ceramic Rotary Valve, Concisely Liquid flow Path
- Various Calibration method, Auto Calibrate
- Traceability Calibration & Control , The results are reliable
- Up to 1000 sample results stored, easy to check

WWW.INSIGNIABIMEDICAL.COM

Specifications

Method	Multi-Biosensor Technology			
Sample Size	60-120 μ l			
Speed	≥ 80 tests/hour			
	Performance	Detection Range(mmol/L)	Precision(CV)	Resolution(mmol/L)
Measure scope & Precision	K ⁺	0.50 ~ 20.00	$\leq 1.0\%$	0.01
	Na ⁺	15.0 ~ 200.0	$\leq 1.0\%$	0.1
	Cl ⁻	15.0 ~ 200.0	$\leq 1.0\%$	0.1
Printer	Inner thermal printer			
Expand Function	RS232 & COM Connector, Computer and Barcode reader available			
Ambient Conditions	0-40 $^{\circ}$ C, Relative humidity $\leq 85\%$			
Power	AC 80V ~ 260V, 50-60Hz, ≤ 80 W			



WWW.INSIGNIABIOMEDICAL.COM

B&E SCIENTIFIC INC.

Add : 804 Thorndale Ave, Bensenville IL 60106 USA

Tel: +1 630 974 5159 Fax: +1 630 974 5160

<http://www.besic.com>